

List of Authorized At-Home OTC COVID-19 Diagnostic Tests

This table lists FDA-authorized at-home OTC COVID-19 diagnostic tests. This table includes links to home use instructions for each test. For more information about each test, including the Letter of Authorization and authorized labeling, see [In Vitro Diagnostics EUAs: Tables of IVD EUAs](#).

Manufacturer and Test Name (Links to Instructions for Use)	Test Type	Who can use this test: Symptoms	Who can use this test: Age	Other Details
Abbott Diagnostics Scarborough, Inc.: BinaxNOW COVID-19 Ag Card Home Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 15 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Requires supervision of a telehealth proctor and a smartphone or computer Results in 15 minutes
Abbott Diagnostics Scarborough, Inc.: BinaxNOW COVID-19 Antigen Self Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 15 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Smartphone optional Results in 15 minutes
Access Bio, Inc.: CareStart COVID-19 Antigen Home Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 14 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Smartphone optional Alternate brand name: On/go COVID-19 Antigen Self-Test Results in 10 minutes
ACON Laboratories, Inc.: Flowflex COVID-19 Antigen Home Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 14 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Smartphone optional Alternate brand name: On/Go One COVID-19 Antigen Home Test Results in 15 minutes
Becton, Dickinson and Company (BD): BD Veritor At-Home COVID-19 Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 14 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Requires a smartphone Results in 15 minutes
Celltrion USA, Inc.: Celltrion DiaTrust COVID-19 Ag Home Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 14 years and older 	<ul style="list-style-type: none"> Mid-turbinate nasal swab Smartphone optional Results in 15 minutes
Cue Health Inc.: Cue COVID-19 Test for Home and Over The Counter (OTC) Use	Molecular	<ul style="list-style-type: none"> People with or without symptoms 	<ul style="list-style-type: none"> Age 18 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Requires a Cue Cartridge Reader (sold separately) Requires a smartphone Results in 20 minutes
Detect, Inc.: Detect Covid-19 Test	Molecular	<ul style="list-style-type: none"> People suspected of COVID-19 People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> Age 14 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Nasal swab Requires a Detect Hub (sold separately) Requires a smartphone Results in 60 minutes
Ellume Limited: Ellume COVID-19 Home Test	Antigen	<ul style="list-style-type: none"> People with or without symptoms 	<ul style="list-style-type: none"> Age 16 years and older Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> Mid-turbinate nasal swab Requires a Smartphone Results in 15 minutes
iHealth Labs, Inc.: iHealth COVID-19 Antigen Rapid Test	Antigen	<ul style="list-style-type: none"> People with symptoms that began within the last 7 days 	<ul style="list-style-type: none"> Age 15 years and older 	<ul style="list-style-type: none"> Nasal swab Smartphone optional

Manufacturer and Test Name (Links to Instructions for Use)	Test Type	Who can use this test: Symptoms	Who can use this test: Age	Other Details
		<ul style="list-style-type: none"> • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Results in 15 minutes
InBios International Inc: SCoV-2 Ag Detect Rapid Self-Test	Antigen	<ul style="list-style-type: none"> • People with symptoms that began within the last 5 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Results in 20 minutes
Lucira Health, Inc.: Lucira CHECK-IT COVID-19 Test Kit	Molecular	<ul style="list-style-type: none"> • People with or without symptoms 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Smartphone optional • Results in 30 minutes
Maxim Biomedical, Inc.: MaximBio ClearDetect COVID-19 Antigen Home Test	Antigen	<ul style="list-style-type: none"> • People with symptoms that began within the last 5 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Results in 15 minutes
OraSure Technologies, Inc.: InteliSwab COVID-19 Rapid Test	Antigen	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 18 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Smartphone optional • Results in 30 minutes
Quidel Corporation: QuickVue At-Home OTC COVID-19 Test	Antigen	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Results in 10 minutes
SD Biosensor, Inc.: COVID-19 At-Home Test	Antigen	<ul style="list-style-type: none"> • People with symptoms that began within the last 6 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Results in 20 minutes
Siemens Healthineer: CLINITEST Rapid COVID-19 Antigen Self-Test	Antigen	<ul style="list-style-type: none"> • People with symptoms that began within the last 7 days • People without symptoms. The test is to be performed two times over three days (serial testing) 	<ul style="list-style-type: none"> • Age 14 years and older • Age 2 years and older when collected by an adult 	<ul style="list-style-type: none"> • Nasal swab • Results in 15 minutes